

Impact Dental Laboratory

200-5300 Canotek Road
Ottawa, Ontario K1J 1A4
www.impact-dental.com

Phone: (613) 746-0602
Toll-free: (800) 668-4691
Fax: (613) 746-5170

Maximum Protrusive Measurement

Enclosed please find a complimentary ruler to conveniently and painlessly take the patient's maximum protrusive measurement.

Using the NTI maximum protrusive ruler, measure from the facial aspect of the Maxillary Centrals to the incisal edges of the Mandibular Centrals.

